



BIO-HAZARD SPILL KIT 25 Ltr

SPLSK25BHZ ABSORBENCY: 25 Ltr

CONTAINS:

> One Set of Perosnal Protective Equipments

> Safety Goggles: 1 Nos

> 3Ply MB Mask: 1 Nos

> Chemical Coverall: 1 Nos

> Safety Gloves: 1 pair

> Baricadding Tape: 1 Nos

> Shoe Cover: 1 Pair

> Clinical Wastage Bag with Tie: 5 Nos

> Absorbent Granules: 2500 Gms

> Caution Mat Yellow: 3 nos

> Chemical Absorbent Pads (Green & Blue): 20Nos

> Chemical Absorbent Socks (8cm x 1.2m): 2 Nos

> Pan & Brush: 1 Nos

SS 304 BOX WITH LUCK: 1 Nos
Rejected X-Ray Sheet: 2 Nos
Disinfectant Liquids: 250 ML

> Blotting Papers: 5 nos

SPILFIX SPLSK100BHZ Biohazard Spill Kits Specially Custome design and packed for quick & Effcient Response against Large Spillage.

- > Put on Protective Gears Supplied with Kit.
- > Sprinkle Granular over the Spill and leave for 2 min.
- > Once Granules have absorbed the spill, collect same using x-ray sheet & Scoop and brush.
- > Sanitise the Spill area with disinfectant liquids.
- > Clean the Spill area using Mats & Pads.
- > Bag and contents must be disposed of in accordance with local waste disposal policies for biohazard waste.
- > As per NABH standard and satisfied ISO guidlines